

---

## Product Datasheet

### PKC theta (Phospho-T538) polyclonal antibody (orb1678124)

**Description**

Liquid in 0.42% Potassium

pH 7.3, 30% glycerol, and 0.01% sodium a...

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

IHC, WB

**Immunogen**

KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PKC theta with a phosphorylation site at T538. The exact sequence is proprietary.

**Preservatives**

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

**Concentration**

1mg/ml

**Storage**

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

**Note**

For research use only

**Clonality**

Polyclonal

**MW**

~ 78 kDa

**Uniprot ID**[Q04759](#)**Dilution Range**

WB (1/500 - 1/1000), IHC (1/50 - 1/200)

**Expiration Date**

12 months from date of receipt.