

## Product Datasheet

### A36R-EEV polyclonal antibody (orb1678396)

<b>Description</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant protein of human A36R-EEV.
<b>Preservatives</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>Q5IXM9</b>
<b>Dilution Range</b>	WB: 1:500~1:2000 ELISA: 1:5000~10000
<b>Expiration Date</b>	12 months from date of receipt.