

Product Datasheet

psmart-1 polyclonal antibody (orb1678405)

Description	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant protein of human psmart-1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q12306
Dilution Range	WB: 1:500~1:2000 ELISA: 1:5000~10000
Expiration Date	12 months from date of receipt.