

## Product Datasheet

# Yellow fever virus E protein ectodomain polyclonal antibody (orb1678424)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2.                  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, WB  |
| <b>Immunogen</b>           | Recombinant protein of YFV-E.  |
| <b>Preservatives</b>       | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.                 |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | AY640589.1   |
| <b>Dilution Range</b>      | WB: 1:500~1:2000 ELISA: 1:5000~10000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |