

Product Datasheet

Recombinant Human AIDS (HIV1) antigen (orb1678434)

Description	Renin is a highly specific endopeptidase that mediates the cleavage of the circulating substrate ang...
Species/Host	Mouse
Tested Applications	CLIA, EIA, LF
Target	HIV
Preservatives	PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4
Concentration	Lot specific(±20%)
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only
Application notes	Method: sandwich method. Use: Coating. Quality control: Mol Mass: 8.36 kDa AP Mol Mass: 8.50 kDa Purity: ≥90% as determined by reducing SDS-PAGE
Isotype	IgG1 kappa
Purity	> 95% as determined by SDS-PAGE.
Source	E.coli
Uniprot ID	P00797
Expiration Date	6 months from date of receipt.