

## Product Datasheet

### AP2S1 Antibody Pair (orb1678556)

<b>Description</b>	AP2S1 Antibody Pair
<b>Reactivity</b>	Bovine, Canine, Equine, Guinea pig, Human, Porcine, Rat, Sheep
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	50% Glycerol, 0.01M PBS, PH 7.4
<b>Form/Appearance</b>	Liquid
<b>Note</b>	For research use only
<b>Application notes</b>	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher.
<b>Isotype</b>	Capture Isotype: IgG. Detection Isotype: IgG
<b>Uniprot ID</b>	<a href="#">P53680</a>
<b>Expiration Date</b>	12 months from date of receipt.