

Product Datasheet

RPS16 Antibody Pair (orb1678564)

Description	RPS16 Antibody Pair
Reactivity	Bovine, Canine, Goat, Guinea pig, Human, Mouse, Rat
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	P62249
Expiration Date	12 months from date of receipt.