

## Product Datasheet

# Human Phosphodiesterase 5A, cGMP Specific (PDE5A) ELISA Kit (orb1678631)

<b>Description</b>	The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of PDE5A in human ser...
<b>Range</b>	1.56-100ng/ml
<b>Storage</b>	Shipped at 2-8°C. Longterm Storage at -20° C
<b>Note</b>	For research use only
<b>Sensitivity</b>	0.59ng/ml
<b>Expiration Date</b>	Please enquire.