

## Product Datasheet

# Human Neuropeptide Y (NPY) ELISA Kit (orb1678646)

<b>Description</b>	The kit is a competitive inhibition enzyme immunoassay technique for the in vitro quantitative measu...
<b>Range</b>	24.69-2000 pg/ml
<b>Storage</b>	Shipped at 2-8°C. Longterm Storage at -20° C
<b>Note</b>	For research use only
<b>Sensitivity</b>	9.44pg/ml
<b>Expiration Date</b>	Please enquire.