

## Product Datasheet

# Human Corticotropin Releasing Hormone (CRH) ELISA Kit (orb1678667)

<b>Description</b>	The kit is a competitive inhibition enzyme immunoassay technique for the in vitro quantitative measu...
<b>Range</b>	12-10,000 pg/ml
<b>Storage</b>	Shipped at 2-8°C. Longterm Storage at -20° C
<b>Note</b>	For research use only
<b>Sensitivity</b>	4.91pg/ml
<b>Expiration Date</b>	Please enquire.