

Product Datasheet

Recombinant Active Cold Inducible RNA Binding Protein (CIRBP) (orb1678756)

Description	Recombinant protein, expressed in Escherichia coli (E. coli) and purified by standard chromatographi...
Preservatives	20mM Tris, 150mM NaCl, pH 8.0, containing 0.01% skl, 5%Trehalose.
Storage	Shipped at 2-8°C. Aliqote and Storet at -80°C; Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Cell culture; Activity Assays. (May be suitable for use in other assays to be determined by the end user.)
Purity	> 95% pure as determined by Silver stained SDS-PAGE analysis.
MW	25KDa
Expiration Date	6 months from date of receipt.