


Product Datasheet

CD315 Antibody (orb1678983)

Description	CD315, also known as prostaglandin F2 receptor negative regulator, is an approximately 135 kDa trans...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	HeLa cells
Target	CD315
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	Q9P2B2
Entrez	5738
Expiration Date	12 months from date of receipt.



Separation of A431 cells stained using a...