

Product Datasheet

Hyponine D (orb1680044)

Description	Hyponine D
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	C ₄₇ H ₅₀ N ₂ O ₁₈
Purity	98%
MW	930.913
CAS Number	259823-31-9
SMILES	<chem>[H][C@@]12[C@@H](OC(=O)c3ccccc3)[C@@]34O[C@@]1(C)COC(=O)c1ccnc1[C@@H](C)[C@H](C)(=O)O[C@@]([H])([C@H](OC(=O)c1ccnc1)[C@H](OC(C)=O)[C@@]3(COC(C)=O)[C@H](OC(C)=O)[C@@H]2OC(C)=O)[C@]4(C)O TLB:15:14:2:50.61.56,16:14:2:50.61.56 </chem>
Expiration Date	12 months from date of receipt.