

## Product Datasheet

### Hyponine E (orb1680045)

<b>Description</b>	Hyponine E
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{45}H_{48}N_2O_{19}$
<b>Purity</b>	98%
<b>MW</b>	920.874
<b>CAS Number</b>	226975-99-1
<b>SMILES</b>	<chem>[H][C@@]12[C@@H](OC(=O)c3ccco3)[C@@]34O[C@@]1(C)COC(=O)c1cccnc1C(C)C(C)C(=O)O[C@@]([H])([C@H](OC(C)=O)[C@H](OC(C)=O)[C@@]3(COC(C)=O)[C@H](OC(C)=O)[C@@H]2OC(=O)c1cccnc1)[C@]4(C)O</chem>   TLB:14:13:2:44.55.50,15:13:2:44.55.50
<b>Expiration Date</b>	12 months from date of receipt.