

Product Datasheet

Peucenidin (orb1680173)

Description	Peucenidin
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{22}O_7$
Purity	98%
MW	386.4
CAS Number	33044-93-8
SMILES	<chem>O(C(C)=O)[C@@H]1C=2C3=C(C=CC2O[C@@H]1C(OC(C=C(C)C)=O)(C)C)C=CC(=O)O3</chem>
Expiration Date	12 months from date of receipt.