

Product Datasheet

Polygalasaponin E (orb1680221)

Description	Polygalasaponin E
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{47}H_{74}O_{19}$
Purity	98%
MW	943.09
CAS Number	882664-72-4
SMILES	<chem>[H][C@@]12CC(C)(C)CC[C@@]1(CC[C@]1(C)C2=CC[C@]2([H])[C@@]3(C)C[C@H](O)[C@H](O)[C@@](C)(C(O)=O)[C@]3([H])CC[C@@]12C)(=O)O[C@@H]1O[C@H](CO)[C@@H](O)[C@H](O)[C@H]1O[C@@H]1O[C@@H](C)[C@H](O[C@@H]2OC[C@@H](O)[C@H](O)[C@H]2O)[C@@H](O)[C@H]1O t:15 </chem>
Expiration Date	12 months from date of receipt.