

Product Datasheet

Isomargaritene (orb1680345)

Description	Isomargaritene
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{28}H_{32}O_{14}$
Purity	98%
MW	592.55
CAS Number	64271-11-0
SMILES	<chem>O([C@H]1[C@@H](O[C@H](CO)[C@@H](O)[C@@H]1O)C=2C(O)=C3C(=CC2O)OC(=CC3=O)C4=CC=C(OC)C=C4)[C@H]5[C@H](O)[C@H](O)[C@@H](O)[C@H](C)O5</chem>
Expiration Date	12 months from date of receipt.