

Product Datasheet

Tagitinin F (orb1681617)

Description	Tagitinin F
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{19}H_{24}O_6$
Purity	98%
MW	348.39
CAS Number	59979-57-6
SMILES	<chem>CC(C)C(=O)O[C@@H]1C[C@@]2(C)O[C@](O)(C=C2)C(C)=C/[C@H]2OC(=O)C(=C)[C@H]12</chem>
Expiration Date	12 months from date of receipt.