

Product Datasheet

Human CST1 Protein (orb168164)

Description	Recombinant of Human CST1 protein
Preservatives	CST1 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 10% glycerol and 100mM NaCl.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMWSPKEE DRIIPGGIYN ADLNDEWVQR ALHFAISEYN KATKDDYRR PLRVLRARQQ TVGGVNYFFD VEVGR TICTK SQPNLDTCAF HEQPELQKKQ LCSFEIYEVW WENRRSLVKS RCQES
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.