

Product Datasheet

beta-Dihydroplumericinic acid (orb1683079)

Description	beta-Dihydroplumericinic acid
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{14}H_{14}O_6$
Purity	98%
MW	278.26
CAS Number	59204-61-4
SMILES	<chem>[H][C@]12O[C@]3([H])[C@H](CC)C(=O)O[C@]33C=C[C@]([H])(C(=CO1)C(O)=O)[C@]23[H] c:13,17 </chem>
Expiration Date	12 months from date of receipt.