

## Product Datasheet

### Humantenine (orb1684539)

<b>Description</b>	Humantenine
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{26}N_2O_3$
<b>Purity</b>	98%
<b>MW</b>	354.44
<b>CAS Number</b>	82375-29-9
<b>SMILES</b>	<chem>[H][C@]12CO[C@@H]3C[C@H]1\C(CN(C)[C@@]2([H])C[C@@]31C(=O)N(OC)c2cccc12)=C\C</chem>
<b>Expiration Date</b>	12 months from date of receipt.