

Product Datasheet

Gymnestrogenin (orb1684554)

Description	Gymnestrogenin
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{30}H_{50}O_5$
Purity	98%
MW	490.72
CAS Number	19942-02-0
SMILES	<chem>C(O)[C@]12[C@](C=3[C@@](C)(C[C@@H]1O)[C@@]4(C)[C@](CC3)([C@]5(C)[C@@](CC4)([C@@](CO)(C)[C@@H](O)CC5)[H])[H])(CC(C)(C)[C@@H](O)C2)[H]</chem>
Expiration Date	12 months from date of receipt.