

Product Datasheet

Ganoderol A (orb1684974)

| | |
|------------------------|---|
| Description | Ganoderol A |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{30}H_{46}O_2$ |
| Purity | 98% |
| MW | 438.69 |
| CAS Number | 104700-97-2 |
| SMILES | <chem>[H][C@@]12CC=C3C(=CC[C@]4(C)[C@H](CC[C@@]34C)[C@H](C)CC\C=C(/C)CO)[C@@]1(C)CCC(=O)C2(C)C</chem> |
| Expiration Date | 12 months from date of receipt. |