

Product Datasheet

APJ receptor agonist 7 (orb1685741)

| | |
|------------------------|--|
| Description | APJ receptor agonist 7 |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{32}H_{46}N_4O_4$ |
| Purity | 98% |
| MW | 550.73 |
| CAS Number | 2100850-28-8 |
| SMILES | <chem>O(C)C1=C(C=2N(N=C(C(N[C@@H](CCCCCCCC3)CC(NC4CCC4)=O)=O)C2)C5CCCC5)C(OC)=CC=C1</chem> |
| Expiration Date | 12 months from date of receipt. |