

Product Datasheet

RORyt agonist 2 (orb1685900)

| | |
|------------------------|--|
| Description | RORyt agonist 2 |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{30}H_{30}F_3N_3O_4S$ |
| Purity | 98% |
| MW | 585.64 |
| CAS Number | 2663787-92-4 |
| SMILES | <chem>C(C)(C)C1N(C=2C(=CC(C(N[C@@H](CC#N)C3=CC=C(S(CC)(=O)=O)C=C3)=O)=CC2)OC1)C4=CC=C(C(F)(F)F)C=C4</chem> |
| Expiration Date | 12 months from date of receipt. |