

Product Datasheet

RORyt agonist 3 (orb1686023)

Description	RORyt agonist 3
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{34}H_{37}N_3O_3S$
Purity	98%
MW	567.74
CAS Number	2664106-24-3
SMILES	<chem>N(CC1CC1)(C2=CC(=CC(=C2)C3CC3)C4CC4)C5=CC=C(C(N[C@@H](CC#N)C6=CC=C(S(CC)(=O)=O)C=C6)=O)C=C5</chem>
Expiration Date	12 months from date of receipt.