

## Product Datasheet

### RORyt inverse agonist 28 (orb1686163)

<b>Description</b>	RORyt inverse agonist 28
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{23}H_{18}Cl_2F_3NO_5S$
<b>Purity</b>	98%
<b>MW</b>	548.36
<b>CAS Number</b>	2741870-21-1
<b>SMILES</b>	<chem>O(C(F)(F)F)C1=C(C=CC=C1)C2=C(Cl)C=C(NC([C@@H](CO)C3=CC=C(S(C)(=O)=O)C=C3)=O)C=C2Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.