

Product Datasheet

LXR agonist 1 (orb1686315)

Description	LXR agonist 1
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{27}H_{26}F_3N_3O_3S$
Purity	98%
MW	529.57
CAS Number	1779524-90-1
SMILES	<chem>C(C)C1=C(C2=C(C(C(F)(F)F)=CC=C2)N=C1)C3=CC(OC4=CC(S(N(CCN)C)(=O)=O)=CC=C4)=CC=C3</chem>
Expiration Date	12 months from date of receipt.