

Product Datasheet

S1p receptor agonist 2 (orb1687527)

Description	S1p receptor agonist 2
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{23}ClN_2O_4$
Purity	98%
MW	438.9
CAS Number	1354908-17-0
SMILES	<chem>C(O)(=O)[C@H]1C[C@@H](N2CC3=C(OC(=N3)C4=CC=C(OCC5=CC=C(C)C=C5)C=C4)CC2)C1</chem>
Expiration Date	12 months from date of receipt.