

Product Datasheet

5-HT2A receptor agonist-2 (orb1688528)

Description	5-HT2A receptor agonist-2
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{20}H_{21}ClN_2O_3$
Purity	98%
MW	372.85
CAS Number	528525-37-3
SMILES	<chem>O(C)C1=C(C2C3=C(C=4C(N3)=CC=C(Cl)C4)CCN2)C=C(OC)C(OC)=C1</chem>
Expiration Date	12 months from date of receipt.