

Product Datasheet

PPAR γ agonist 3 (orb1688570)

| | |
|------------------------|--|
| Description | PPAR γ agonist 3 |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | C ₂₄ H ₂₃ N ₃ O |
| Purity | 98% |
| MW | 369.46 |
| CAS Number | 2011801-48-0 |
| SMILES | <chem>O=C(N)C=1C=CC=CC1C=2C=CC(=CC2)CN3C(=NC=4C=CC=CC43)CC</chem> C |
| Expiration Date | 12 months from date of receipt. |