

Product Datasheet

Kizuta saponin K11 (orb1690110)

Description	Kizuta saponin K11
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{61}H_{98}O_{27}$
Purity	98%
MW	1263.428
CAS Number	97240-03-4
SMILES	<chem>[H][C@@]12CC(C)(C)CC[C@@]1(CC[C@]1(C)C2=CC[C@]2([H])[C@@]3(C)CC[C@H](O)[C@@H]4OC[C@H](O)[C@H](O)[C@H]4O[C@@H]4O[C@@H](C)[C@H](O)[C@@H](O)[C@H]4O)[C@@](C)(CO)[C@]3([H])CC[C@@]12C)C(=O)O[C@@H]1O[C@H](CO[C@@H]2O[C@H](COC(C)=O)[C@@H](O[C@@H]3O[C@@H](C)[C@H](O)[C@@H](O)[C@H]3O)[C@H](O)[C@H]2O)[C@@H](O)[C@H](O)[C@H]1O t:15 </chem>
Expiration Date	12 months from date of receipt.