

Product Datasheet

APJ receptor agonist 1 (orb1691035)

Description	APJ receptor agonist 1
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{31}H_{26}ClN_3O_3$
Purity	98%
MW	524.02
CAS Number	2287153-38-0
SMILES	<chem>CCCC(NC(=O)c1ccc(c(c1)C(O)=O)-c1c(Cl)cccc1-c1nc2cccc2[nH]1)c1cccc1</chem>
Expiration Date	12 months from date of receipt.