

Product Datasheet

C3a Receptor Agonist (orb1692142)

| | |
|------------------------|---|
| Description | C3a Receptor Agonist |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{27}H_{35}N_3O_2$ |
| Purity | 97.50% |
| MW | 433.59 |
| CAS Number | 944997-60-8 |
| SMILES | <chem>C(C(NC1CCN(C(CCC=2C=CC=NC2)=O)CC1)=O)(C3CCCCC3)C4=CC=CC=C4</chem> |
| Expiration Date | 12 months from date of receipt. |