

Product Datasheet

Glucocorticoid receptor agonist-1 (orb1692476)

| | |
|------------------------|--|
| Description | Glucocorticoid receptor agonist-1 |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{35}H_{39}NO_6$ |
| Purity | 98% |
| MW | 569.69 |
| CAS Number | 2166375-82-0 |
| SMILES | <chem>C(CO)(=O)[C@]12[C@]3(C)[C@@](C[C@]1(O[C@H](O2)C4=CC=C(CC5=CC(N)=CC=C5)C=C4)[H])([C@]6([C@]([C@@H](O)C3)([C@]7(C)C(CC6)=CC(=O)C=C7)[H])[H])[H]</chem> |
| Expiration Date | 12 months from date of receipt. |