

Product Datasheet

Dengue 1 NS1 antibody (orb169705)

Description	Rabbit polyclonal antibody to Dengue 1 NS1
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	0.2M glycine pH 7.5 and 0.02% NaN ₃ .
Storage	Stability: Should be stored below -18°C. Please prevent freeze thaw cycles
Note	For research use only
Application notes	Applications: Immunoassay
Clonality	Polyclonal
Purity	Protein is > 95% pure as determined by 12% PAGE (coomassie staining). Purified by protein A chromatography.
Source	Rabbit
Expiration Date	12 months from date of receipt.