

Product Datasheet

Dengue 2 NS1 antibody (orb169706)

Description	Rabbit polyclonal antibody to Dengue 2 NS1
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	PBS + 0.02% sodium azide.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Polyclonal Antibodies, Anti-Viral
Clonality	Polyclonal
Purity	Protein is > 90% pure.
Expiration Date	12 months from date of receipt.