

## Product Datasheet

### Biperiden lactate (orb1698140)

<b>Description</b>	Biperiden lactate
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{24}H_{35}NO_4$
<b>Purity</b>	98%
<b>MW</b>	401.54
<b>CAS Number</b>	7085-45-2
<b>SMILES</b>	<chem>O=C(O)C(O)C.OC(C=1C=CC=CC1)(CCN2CCCCC2)C3CC4C=CC3C4</chem>
<b>Expiration Date</b>	12 months from date of receipt.