


Product Datasheet

Mouse CCL28 ELISA Kit (orb170143)

Description	Mouse CCL28 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: I23-R130
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 400 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	14570 MW
Uniprot ID	Q9JIL2
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CCL28
PicoKine
ELISA Kit
standard ...