


Product Datasheet

Rat FGF1/Fgf Acidic ELISA Kit (orb170155)

| | |
|----------------------------|---|
| Description | Rat FGF1/Fgf Acidic ELISA Kit |
| Reactivity | Rat |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: F16-D155 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:10, concentration in serum and plasma is around 5.6 ng/ml |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| MW | 17418 MW |
| Uniprot ID | P61149 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Rat FGF1 PicoKine ELISA Kit standard cur...