


Product Datasheet

Mouse Progranulin ELISA Kit (orb170168)

Description	Mouse Progranulin ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T18-L589
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is 492-2197 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	63458 MW
Uniprot ID	P28798
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse Progranulin PicoKine ELISA Kit sta...