


## Product Datasheet

### Mouse Betacellulin/BTC ELISA Kit (orb170169)

<b>Description</b>	Mouse Betacellulin/BTC ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	7.8 pg/ml - 500 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: D32-Y111
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 7.8 pg/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	19664 MW
<b>Uniprot ID</b>	<b>Q05928</b>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse Betacellulin/BTC PicoKine ELISA Kit...