


Product Datasheet

Human Cystatin B/CSTB ELISA Kit (orb170179)

Description	Human Cystatin B/CSTB ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M2-F98
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 1.7 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	11140 MW
Uniprot ID	P04080
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Cystatin B PicoKine ELISA Kit stan...