


Product Datasheet

Human PSP94/MSMB ELISA Kit (orb170181)

Description	Human PSP94/MSMB ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S21-I114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	12865 MW
Uniprot ID	P08118
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
PSP94
PicoKine
ELISA Kit
standard ...