

Product Datasheet

Fexarene (orb1702272)

Description	Fexarene
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{32}H_{33}NO_3$
Purity	98%
MW	479.62
CAS Number	574013-68-6
SMILES	<chem>COC(=O)C=Cc1cccc(c1)N(Cc1ccc(C=Cc2ccccc2)cc1)C(=O)C1CCCCC1</chem> w:4.3,18.18
Expiration Date	12 months from date of receipt.