

## Product Datasheet

### ZDHHC17 Immunizing peptide (orb17027)

<b>Description</b>	ZDHHC17 Immunizing peptide
<b>Preservatives</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Storage</b>	Shipped at ambient temperature, store at - 20°C
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb17027 (ZDHHC17 antibody)
<b>Protein Sequence</b>	CYDQISGSGYQLV
<b>Expiration Date</b>	12 months from date of receipt.