

Product Datasheet

Tankyrase Inhibitors (TNKS) 49 (orb1704012)

Description	Tankyrase Inhibitors (TNKS) 49
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{24}N_4O_3S$
Purity	98%
MW	448.54
SMILES	<chem>O=C1NC(SCCC(N[C@H]2CC[C@H](OC3=CC=C(C#N)C=C3)CC2)=O)=N C4=CC=CC=C41</chem>
Expiration Date	12 months from date of receipt.