

Product Datasheet

RORyt Inverse agonist 2 (orb1708650)

| | |
|------------------------|--|
| Description | RORyt Inverse agonist 2 |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{27}H_{25}F_8NO_5S$ |
| Purity | 98% |
| MW | 627.54 |
| CAS Number | 1801197-92-1 |
| SMILES | <chem>S(=O)(=O)([C@@]1(CN(C(=O)[C@H]2CC[C@H](C(O)=O)CC2)CC1)C3=C=C(C(C(F)(F)F)(C(F)(F)F)F)C=C3)C4=CC=C(F)C=C4</chem> |
| Expiration Date | 12 months from date of receipt. |