

Product Datasheet

RORyt Inverse agonist 3 (orb1708652)

Description	RORyt Inverse agonist 3
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{29}H_{31}Cl_2N_5O_4S$
Purity	98%
MW	616.56
CAS Number	2364429-77-4
SMILES	<chem>O(CC1=C(Cl)C(S(=O))(=O)N2[C@H](C(NC)=O)CCC2)=CC=C1Cl)C=3C4=C(C(=CC(C(C)C)=N4)C=5N(C)C=NC5)C=CC3</chem>
Expiration Date	12 months from date of receipt.