

Product Datasheet

GLP-1 receptor agonist 3 (orb1709603)

Description	GLP-1 receptor agonist 3
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{31}H_{30}FN_5O_4$
Purity	98%
MW	555.6
CAS Number	2230200-09-4
SMILES	<chem>C(N1C=2C(N=C1CN3CCC(CC3)C=4N=C(OCC5=C(F)C=C(C#N)C=C5)C=CC4)=CC=C(C(O)=O)C2)[C@H]6CCO6</chem>
Expiration Date	12 months from date of receipt.